



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

## AIR CURTAIN

- Wooden with sun-mica/MS sheet duly powder coated
- Size Available 3' feet / 4' feet / 5' feet
- Powder consumption 220/230 V. AC supply single phase, 50 Hz.

**IIC-150**




## BOTTLE WASHING MACHINE

- Mounted on heavy duty angle iron stand duly painted with tank made of G.I. / S.S.
- Electrically motor operated, with quick released chuck which can accommodate of different variety of nylon brushes for cleaning test tube bottles etc., with 12 rinses to drill chuck.

### Available In Size

- |  |            |
|--|------------|
| With G.I. soaking tank of size 30"x20"x12" - | IIC 151 A. |
| With G.I. soaking tank of size 30"x20"x18" - | IIC 151 B. |
| With S.S. soaking tank of size 30"x20"x12" - | IIC 151 C. |
| With S.S. soaking tank of size 30"x20"x18" - | IIC 151 D. |

**IIC-151**



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in