



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)



BIO SAFETY CABINET

Available in
Class-II, Type-A2 .
Class-II, Type B2 (Total Exhaust Model)

Bio Safety Cabinet Class-II

Working Table of S.S. 304 grade. Body of construction can be of Wodden/Mild Steel/Stainlesss Steel

- HEPA of 0.3 u of an efficiency of 99.97% down to 0.3u media is of glass fiber & made of pleated back & frow made of aluminium foil & packed in a anodized aluminium frame.
- Pre filter of dry fiber washable type with the frame.
- Air flow: 70% recirculation & 30% exhaust.
- Statically balanced motor blower assembly
- Velocity at output of HEPA is 90+20 FPM
- Static Pressure Manometer
- Transparent front Door.
- Normal working fluorescent light
- UV light 1X 15/30 watts.
- Gas cock

Options:-

- (A) Virus Burnout System
- (B) HEPA filter in Exhaust
- (C) Magnehelic Gauge
- (D) Glove Port

Size:- 2'X2'X2'

3'X2'X2'

4'X2'X2'

6'X2'X2'

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in