



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

**IIC-140**



## DEEP FREEZER / LOW TEMP. VERTICAL/HORIZONTL MODEL (ECO FRIENDLY)

- Refrigerated Low Temperature cabinet.
- Better puff Insulated to avoid cooling loss.
- Silent energy saving compressor.
- Polish 304 grade Stainless Steel sheet interior.
- Outer MS/ CRC sheet powder Coated.
- Optional : Outer S. Steel Body against customer requirement
- Provide with Digital Temp. Controller Cum Indicator,/Thermostat Bi-Metallic.
- Door Front/ Top opening for Horizontal model duly Puf insulated.

### Applications :

Agriculture, Bio Tech, Rubber, Life Science , Pharmaceutical, Research Lab. Quality inspection.

Model	IIC 140	IIC 140 A	IIC 140B
Capacity (L)	113 Liter	173 Liter	282 Liter
Interior mm	455 x 610 x 410	510 x 830 x 410	570 x 870 x 570
Shelves	2	2	2/3

### Technical Parameters :

- Temperature Range : Up to - 20°C / 30°C / 40°C
- Temperature Stability : ± 1°C
- Material : CRC Sheet duly powder coated or Stainlesssteell
- Power consumption : 220/230V AC Single Phase 50HZ
- Compressor : Kirloskar /Tecumseh/Imported as per Model requirement
- Refrigerated : R 22/ R 134
- Duty Cycle : Continuous
- Interior : S. Steel 304 Grade polish sheet
- Optional : Voltage Stabilizer ( On Demand)
- Caution : -40 should install and run in A.C. Room

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in