



ICON INSTRUMENTS COMPANY
(Manufacturer of Laboratory Equipment)



FLAME PHOTOMETER

Model No. IIC 414

2 Filter Na, K

Range from 0 to 100 ppm

Narrow band interference glass filter

LPG & dry oil free air

Measure one element at a time

Single Display 2 1/2 digit bright red led Display

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in