



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-175

GLASS WARE DRIER

Dry laboratory glassware quickly and safely with this bench or wall mounted storage unit. By providing continuous warm air circulation through hollow n drying pegs, lab-ware is never too hot to handle. Three sizes of drying pegs accept glassware with large or small bores, and the polyurethane foam filter keeps warm air dust free, for rapid & convenient drying of Laboratory glassware by blower forced air through 12 brass nickel plated pipes fitted with heating elements blower through powerful motor On/Off switch to work on Single phase 220/230 volt A.C.



IIC-176

KJELDAHL DIGESTION & DISTILLATION UNIT

These Kjeldahl digestion & distillation units are very accurate and supreme in functioning for the decomposition of the organic materials and are also used as distillation units. These units are provided a very accurate and robust construction which makes them more stable and reliable for its use. We are Offering these Kjeldahl digestion & distillation units to our clients at a very comfortable range of prices.

Model No.	IIC-176	IIC-176 A
Capacity	Unit of 03 Test	Unit of 06 Test
Heating System	Mantle Type	Mantle Type



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in