



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-116



HEATING MANTEL

- Small & Large capacity Mantle.
- Electrical Heated hard knitted glass yarn without joint.
- One pcs Spun body made from non-rusting Alum. Duly powder coated & stove enamel painted.
- Full feature panel with ER / Temp. Controller, Temp. Indicator & Mains Switch.
- Better mineral / Glass-wool Insulated.
- Useful for heating inflammable liquids.

Applications :

For Research Lab. Pharmaceutical Chemical /Glass Plants, Bio Tech Lab. Educational Institute.

Model	IIC 116	IIC 116 A	IIC 116 B
Capacity (L)	250ml, 500ml	1 Liter, 2 Liter	3, 4, 5 Liter
Heater Wattage "W"	200W, 200W	300W, 450W	2 x 300W 2 x 500W 2 x 500W

Technical Parameters :

- Temperature Range : Ambient to 350°C
- Power Consumption : 220/230V AC ± 10%, 50HZ
- Material : One Pcs Alum. Spun body duly painted / Powder Coated
- Insulation : Mineral / Glass-wool
- Optional : Double Colour body (On Customer Request)
- Control Panel : ER, Temp. Indicator & Mains Switch

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • **E-mail :** iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in