



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

**IIC-126**



## LABARTORY HOT PLATE

- Unique design of outer MS body having powder coated with gap of cast Iron top from the body.
- Kanthal A-1 superior quality coil shaped heater duly embedded in the cast Iron plate to avoid. heat loss having Cera Blanket Insulation.
- All electrical connection are fire proof having sleeve for Heat/Temperature Protection
- Temp. Range is from Ambient to 300°C.
- Temp. Controlled by thermostat controller
- Cast Iron top is having 7.5" Diameter having heat resistant Paint.
- Optional : Can be supplied with Digital controller in place of Thermostat.

### Applications :

Laboratory use : for educational Institute Hospital, Research Lab. Bio Tech Lab. & Pharmaceutical Industries

Model	IIC 126	IIC 126 A	IIC 126 B	IIC 126 C	IIC 126 D
Size	7.5" Dia	10" x 12" x 6" Rectangular	16" x 10" x 6" Rectangular	18" x 12" x 6" Rectangular	24" x 18" x 6" Rectangular
Heater "W" wattage	1.5 KW.	1.5 KW.	1.5 KW.	2.0 KW.	3.0 KW.

Available model : with ER & with Digital Controller

### Technical Parameters :

Control Panel	:	Energy Regulator, Temp. Indicator & Load Indicator
Material	:	CRC/MS Sheet duly Powder Coated
Optional	:	All S. Steel on request, Digital in place of Thermostate
Duty Cycle	:	Continuous
Power Consumption	:	220/230V AC Supply 50HZ single Phase

**Works :** Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

**Phone :** +91-11-26383676 • **E-mail :** iconinstruments.india@gmail.com, sales@iconinstruments.co.in

**Website:** www.iconinstruments.co.in