



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

HUMIDIFIER

- It is of splash plate rotating disc type.
- Designed to give as high as 95% RH in enclosed area.
- Capable of Atomising up to 3.0 Litre of water/Hr at full evaporation in the driest Ambient Condition.
- Body & Tank made of S. Steel.

Applications :

Pharmaceutical, Quality Lab., Institute

Model	IIC 131	IIC 131 A
Capacity (Liter)	3 Liter	5 Liter

Technical Parameters :

- Material : All S. Steel with castor wheel
- Power Consumption : 220/230V AC Supply 50HZ single Phase
- Duty Cycle : Continuous Interior / Exterior All S. Steel 304 Grade
- Optional : Humidity Controller with Box for above unit o Timer Box for Above unit




TLC KIT (THIN LAYER CHROMATOGRAPHY KIT)

- Alum/ Brass Applicator, As per requirement.
- Thin layer drawing rac.
- Developing chamber.
- Aluminium carrier, spotting Jug, Silica Gel.
- Pipette 0.1 ml.
- Square / Rectangular glass plate & glass spray man.

Applications :

For Research Institute Lab. Medical Institute Collages.



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in