




# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

## IIC-163

### JUICE PRESS (TO EXTRACT JUICE FROM CRUSHED FRUITS)

- Hand-operated unit available with or without a hydraulic pump of capacity 5 tones.
- Mounted on very heavy base and pillars along with hand screw arrangement.
- Supplied with ETEAKWOOD circular plank held together with stainless steel strips.
- Thick M.S. base tray lined with stainless steel is provided with the unit.
- Press for extracting juice from apples, pineapples, carrots, ginger, grapes etc.
  1. Hand screw bucket press
  2. Hydraulic Juice Press.



## IIC-164

### STEM JACKETED KETTLE

- Used for batch heating and cooking of various food products.
- Mounted on a heavy duty mild steel stand with tilting arrangement.
- Jacket and pan both are made of high quality S.S.-304 grade stainless steel of heavy gauge.
- Ten pan is 2/3 jacket for maximum stem utilization and efficiency.

#### Available In Size

- 10 gallons (45ltrs.) - IIC 164 A.
- 25 gallons (110 ltrs.) - IIC 164 B.
- 50 gallons (225 ltrs.) - IIC 164 C.



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in