



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-124



LAMINAR AIR FLOW / HORIZONTAL & VERTICAL CLEAN BENCH

- Better Sunmica & Plywood duly varnished duly painted Excellent wooden & Castor wheel for easy movement.
- Silent Air Blower unique design of Circulation for uniform air movement.
- A Horizontal uni flow clean bench has an open work table featuring a reliable & easy control.
- Low noise & low vibrations motor is adopted for, the speed can be adjusted step by step, The thin switch on panel make a reliable & easy control.
- The work table pattern is made of S. Steel sheet nice & durable. The external panel is finished with sun-mica sheet nicely fixed.
- Full feature Hepa Filters, Pre-filters comprising of non-woven micro fibre glass pleated media, pleats are separated by corrugated Alum. Foil separators, HEPA filters are protect by plastic grill.
- The side panels are made of good quality Thick Perspex sheet duly framed in anodized Alum.
- Working Bench is fitted with Fluorescent Light & UV Germicidal light having separate control in front.
- Voltage required 220/230V AC Supply 50 HZ

Applications :

Laminar Air Flow / Horizontal & Vertical Clean Bench : IIC-124

Outer is finished with sunmica & inner epoxy painted, Working table made of S.S.

Model	IIC 124	IIC 124 A	IIC 124 B	IIC 124 C
Working Zone	2' x 2' x 2'	3' x 2' x 2'	4' x 2' x 2'	6' x 2' x 2'

Outer is made of Cold Rolled steel duly powder coated, Working table made of S.S.

Model	IIC 124-1	IIC 124 -1A	IIC 124 - 1B	IIC 124 - 1C
Working Zone	2' x 2' x 2'	3' x 2' x 2'	4' x 2' x 2'	6' x 2' x 2'

Outer is made of Stainless Steel, Working table made of S.S.

Model	IIC 124-2	IIC 124 -2A	IIC 124 - 2B	IIC 124 - 2C
Working Zone	2' x 2' x 2'	3' x 2' x 2'	4' x 2' x 2'	6' x 2' x 2'

Vertical Model Laminar Air Flow Bench : All above size available A Vertical uniflow clean Bench has a close worktable, Preventing the Air from coming in all other specification i.e. size, No. of Hepa Filters, UV light, Pre Filter motor & Fluorescent Tube remains same as in horizontal model.

- Optional :**
- Magnehelic Gauge to track filter pressure
 - Sun-mica Top instead steel of S. S. Top.

TOTAL EXHAUST MODEL ALSO AVAILABLE WITH SPECIAL ARRANGEMENT FOR NEAR 100% EXHAUST OF AIR.

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in