



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

**IIC-139**



## MELTING POINT

- Useful for determination of Single & mixed melting point.
- Temperature Range is upto 360°C.
- Feature low/high Provision Energy Regulator.
- Aluminum Casted body.
- Outer casing MS sheet duly painted having provision for capillary & Thermometer for testing.
- Supplied with Glass capillary.
- Also Available in Digital Model

## Applications :

For PVC, Rubber Educational, Research Institute

Model	IIC 139
Interior	Alum. Block
Wattage "W"	10 Wattage Appox.

## Technical Parameters :

Temperature Range	: Ambient to 360°C Temp. Control ER (By energy Regulator) Interior Alum. Block duly machined
Material	: CRC/MS sheet duly painted
Duty Cycle	: Continuous
Power Consumption	: 220/230V AC Single Phase 50HZ
Control Panel	: ER Mains & Heating Indicator, Mains

**Works :** Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

**Phone :** +91-11-26383676 • **E-mail :** iconinstruments.india@gmail.com, sales@iconinstruments.co.in

**Website:** www.iconinstruments.co.in