



ICON INSTRUMENTS COMPANY

IMPACT STUDENT MICROSCOPE (IRIS DIAPHRAME TYPE)

with triple revolving nose piece, horse shoe, case separate coarse an fine motion knobs having stopper for slides safety, sub stage condenser N.A. 1.2 with iris diaphragm and filter holder & reflecting mirror alongwith Achromatic Objectives 10X & 45X ,Huygenian Eye Piece 10X & 15X



IMPACT MEDICAL RESEARCH MICROSCOPE

triple revolving nose piece, separate coarse and highly sensitive graduated fine focusing knobs, slide safety stopper, sub stage condenser N.A 1.2 moving up & down through rack & pinion, fixed stage size 120x120mm with precision mechanicals stage & reflecting Mirror ,alongwith Achromatic Objectives 10X ,45X 100X oil , Huygenian Eye Piece 10X & 15X



IMPACT MONOCULAR RESEARCH MICROSCOPE

nose piece, separate coarse and highly sensitivity graduated fine focusing knobs, observation tube inclined at 45 degree and rotatable in 360 degree, fixed stage size 120x140mm with precision mechanicals stage for X & Y movement sub stage bright field condenser N.A 1.2 which is centrable in mount, Built in base illumination system having 6V - 20W HALOGEN LAMP alongwith Achromatic Objectives 10X ,45X 100X oil ,Huygenian Eye Piece 10X & 15X



IMPACT BINOCULAR RESEARCH MICROSCOPE

A latest and improved wide field Binocular Research Microscope with binocular head to ensure highly accurate collimation from stain free observation, Adjustable interpupillary and dioptic distance, heavy duty quadruple revolving nose piece fixed co-axial mechanical stage both sides graduated for slide manipulation, separate coarse and highly sensitive graduated fine focusing knobs sub stage bright field condenser N.A 1.2 which is centrable in mount, built in base HALOGEN illumination system 6V-20W alongwith Achromatic Objectives , 5X ,10X ,45X (SL) & 100X oil , Huygenian Eye Piece 10X & 15X



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in