



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

**IIC-112**

## ROTARY SHAKER

- Better-fabricated assembly frame heavy mild steel section duly painted or powder coated.
- Counter Balanced drive mechanism ensure high stability & reliability in regular operation
- Available Fixed & Variable Speed.
- Maintenance free Induction drive motor.
- Variable speed from 100 to 180 RPM & 20 to 180 RPM as per requirement.
- Uniform platform to accommodate interchangeable clamps of different size for different capacity of Flask.
- Uniform circulation of 19mm diameter to Load platform in Horizontal Plane.
- Uniform Circulation motion of 25mm dia in horizontal plane for heavy duty Rotary Shaker.



## Applications :

Useful mixing /developing of culture, solvent Chemicals. Pathological & various application.

Model	IIC 112	IIC 112 A	IIC 112 B	IIC 112 C
Dimension	12" x 12" Variable Speed	18" x 18" Variable Speed	25" x 25" with Variable Speed 1/4HP Motor	35" x 35" with Fixed Speed 1/2HP Motor
Capacity	9 Flask of 250ml	16 Flask of 250ml	25 x 250ml or 36 x 100ml	36 Flask of 250ml
No. of Tray	01 No.	01 No.	01 No.	01 No.

## Technical Parameters :

- Material : Fabricated assembly from heavy mild steel (MS) Section duly paint/ powder Coated.
- Duty Cycle : Continuous
- Optional : o Digital RPM meter with DC drive @ 5,000
- Power Consumption : 220/230V AC supply ± 10%, single Phase 50HZ

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in