



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-114



SHAKING MACHINE / BOTTLE SHAKER (RECIPROCATING KHAN'S TYPE)

- CRC /MS sheet housing duly painted.
- Fitted with 1/4HP heavy duty motor.
- No. of oscillation fixed at 200 ± 5 per minute.
- Supplied with Top Platform rod on/off switch.
- With carrier for flask, Test tube Rack bottles etc.
- Speed of motor unaffected by normal line voltage fluctuations or loading.

Applications :

For Hospital, Research Lab. Pharmaceutical Institute

Model	IIC 114		IIC 114 A	
Capacity mm	20" x 13" x 4" (500 x 325 x 100mm)		31" x 17" x 4" (775 x 425 x 100mm)	
Motor (H. P.)	1/4 HP		1/4 HP	
Type	Open Body	Cover Body	Open Body	Cover Body

Technical Parameters :

- Material : CRC / MS angle sheet duly painted
- Power Consumption : 220/230V AC Supply 50 HZ Single Phase
- Duty Cycle : Continuous
- Optional : Timer can be fitted (On request)

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • **E-mail :** iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in