



ICON INSTRUMENTS COMPANY

SHERRINGTON'S ROTATING DRUM

- Kymograph: Kymo means wave; Graphic means write. Kymograph is thus the name given to any instrument which records movement on moving surface. The degree of movement can be analyzed from the height of the wave and its relation to time is measured by horizontal distance. It consists of following parts: a) Electric Motor: It runs on 220 volts AC b) Shaft: It is connected to the motor vertically and can be rotated at different speeds.

IIC-330



HAND GRINDING MILL (TABLE CLAMPING)

- Table fixing with hardened grinding discs, provision for course and fine grinding

IIC-331



ASEPTIC CABINET

- Size 95 x 50 x 50 cms. M.S. structure with slopping front, acrylic front and side doors (fully visible, transparent), 3' germicidal ultra-violet tube and 2' fluorescent tubes for white light, with necessary choke blast, starter and switches. Also mains lead; plug and two holes with rings and sleeves for insertion of arms are provided

IIC-332



CONICAL PERCOLATOR

- Made of Stainless Steel capacity 1000ml with lid, drain cock, jointless stainless steel sieve on 3 legged stand

IIC-333



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in