



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

## IIC-122

### SOXHLET EXTRACTION APPTS.

- Small & Large size extraction.
- Electrical heating hard knitted glass yarn.
- MS/CRC sheet outer body having provision for MS top duly painted for 3 & 6 Mantles.
- Supplied with 2 Nos. Vertical / Horizontal rods with gutkas & screw to hold flask.
- Suitable to heat Soxhlet flask of 50 to 500ml capacity.
- Supplied without Glass Parts.



### Applications :

For Research Lab. Pharmaceutical Chemical / Glass Paints educational Institute

Model	IIC 122	IIC 122 A
Capacity (L)	Unit 3 Sets	Unit 6 Sets
Heater "W" Wattage	200 /250W	200 /250W

### Technical Parameters :

- Temperature Range : Ambient to 300°C
- Power Consumption : 220/230V AC supply 50 Hz
- Material : MS/CRC sheet duly powder coated
- Insulation : Mineral / Glass-wool
- Optional : Double Colour body / Top (On Customer request)
- Control Panel : ER 3/6 Nos. (As per requirement) 6 Nos. Toggle Switch 6 Nos. Temperature Indicator

## IIC-123

### KJELDHAL DIGESTION UNIT

Mental type (Without glass parts & Clamps) Useful for estimation of nitrogen contents, The housing of mild steel sheet duly powder coated unit of 3/6 Test as per detail (Given Above with lead fume duct & pair of hangers. Complete in all respect  
 Unit of 3 Sets  
 Unit of 6 Sets



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in