



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-158

LABORATORY STIRRER

- The motor stirrer having maximum speed up to 4000 r.p.m unit is clamped to a C.P. rod fitted on a heavy cast iron base.
- The height and angle of the stirrer can be adjusted according to the requirement.
- Variable speed Motor with speed regulator is fitted with S.S. Stirrer rod to give smooth and noiseless operation

Available In Size

Without Drill Chuck - IIC 158.

With Drill Chuck - IIC 158 A

Optional

Automatic 0-60 minute "GERMAN" Timer.

Digital Model also Available



IIC-159

EMULSIFIER (HOMOGENIZER)

- Fitted with F.H.P. motor having speed regulator and Maximum speed 4000 RPM with perforated S.S. disc.
- The unit is having arrangement for moving in horizontal and vertical direction (Complete on stand).

Available In Size

With 1/20 H.P. Motor and S.S. Disc and speed regulator - IIC159 A.

-do- with 1/4 H.P. Heavy Duty Motor and speed regulator - IIC 159 B.

-do- with 1/2 H.P. Heavy Duty Motor and speed regulator - IIC 159 C.



Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in