



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-152

TEST SIEVES (COARSE SERIES)

- Made of perforated steel plate with square holes.
- Available in following aperture of 125, 160, 100, 90, 80, 75, 63, 53, 50, 45, 40, 37.5, 31.5, 26.5, 25, 22.4, 19, 16, 13.2, 12.2, 11.2, 10.2, 10, 9, 5, 8, 6.7, 6.3, 5.6, 4.75, 2.8 mm



Available In Size

- With brass frame 8" dia.
- With G.I. frame 12" dia.
- With G.I. frame 18" dia.

IIC-153

TEST SIEVES (FINE SERIES IN 8" DIA. BRASS FRAME)

B.S.S.	A.S.T.M.	I.S.I.	Microns	mm.
4	5	480	----	4.75
5	6	340	1812	3.35
6	7	280	2411	2.80
7	8	240	2207	2.36
8	10	200	2057	2.00
10	12	170	1698	1.70
12	14	140	1005	1.40
14	16	120	1200	1.18
16	18	100	1000	1.00
18	20	85	850	0.850
22	25	70	710	0.710
25	30	60	600	0.600
30	35	50	500	0.500
36	40	40	425	0.425
44	45	35	355	0.355
52	50	30	300	0.300
60	60	25	250	0.250
72	70	20	212	0.212
85	80	18	180	0.180
100	100	15	150	0.150
120	120	12	125	0.125
150	140	10	106	0.106
170	170	9	90	0.90
200	200	8	75	0.75
240	230	7	63	0.63
300	270	6	53	0.53
325	300	5	45	0.45
350	325	4	44	0.044
400	380	3	40	0.040
500	---	---	25	0.025



Receiving Pan and Lid made of Brass 8"

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in