



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

**IIC-172**



## TEST TUBE WARMER

- Test tube Warmer has been developed to provide low dry heat. The unit maintains accurately the predetermined temperature. The value of set temperature is totally user settable & is indicated digitally on the front panel of the temperature controller. It has got fully SS body, temperature accuracy is +/- 0.2 Deg C with three anodized aluminum blocks.
- Guidelines & Operating Instruction:
  - 1) Complete unit consist of : S.S. Body, Three Aluminum Blocks, One Temperature Controller & One Power cord.
  - 2) Connect the warmer to earthed power supply, 230VAC 50 Hz using power cord provided.
  - 3) Switch the mains supply in ON position, display will illuminate immediately.
  - 4) Place the blocks on the heated portion.
  - 5) In half an hour, block will attain the desired set temperature (37°C)
  - 6) Warmer is easy to use.

**Works :** Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

**Phone :** +91-11-26383676 • **E-mail :** iconinstruments.india@gmail.com, sales@iconinstruments.co.in

**Website:** www.iconinstruments.co.in