



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

IIC-128

THERMAL STABILITY TEST APPTS.

- Double walled useful for cable, PVC Rubber Testing,
- Solid Alum. Round shape block having no. of holes for Thermal tubes.
- Strip type deluxe quality Kanthal heater to give uniform Temperature.
- S. S. 304 grade S. Steel top having provision for holes for different tubes.
- Full feature Microprocessor based PID Temp. Controller cum indicator to give better accuracy.
- Better Cera blanket Insulated to give minimum heat loss. Avoid heat loss & gives better results.
- Unique handle arrangement to carry equipment during / After Test.



Applications :

Useful for PVC cable & Rubber Industries Research Lab.

Model	IIC 128
Interior	Alum. Block heavy duty
Heater "W" Wattage	350W

Technical Parameters :

Temperature Range	: Ambient to 210°C, working at 200°C continuous
Temperature Accuracy	: $\pm 1^\circ\text{C}$ or better
Digital Setting	: 0.1°C
Digital Resolution	: 0.1°C
Material	: CRC/MS sheet duly powder coated
Interior	: Alum. Block
Duty Cycle	: Continuous
Power Consumption	: 220/230V AC Supply $\pm 10\%$ 50 HZ

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in