



# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

**IIC-174**

## TISSUE CULTURE RACK

- Tissue culture racks are very useful for various types of tissue culture plants & other tissue culture work in tissue culture room.
- Tissue racks are made of tubular mild steel pipes duly painted for longer life.
- Each rack of size 48"x12" (depth) is covered with thick unbreakable plexi glass/ sun mica paste duo board. Which are termite and insect proof fire and weather resistant. Distance between each rack is 15"
- Four fluorescent tubes are properly illuminated by diffused; glare free fluorescent lights are built in upper portion in worktable.
- Are fitted inside the rack/ work area.
- All fluorescent tubes chocks & ultra violet germicidal tube chocks are fitted outside of rack. If desired sodas to avoid heat dissipation inside the rack.
- 0-24 Hours timers are fixed front of tissue culture rack.
- All control panels are fitted tops of tissue culture racks. ON/OFF switch of fluorescent tubes & UV tubes. All control panels are fitted separately in each rack with indicator light.
- The whole rack is provided with four nos. castor wheels at the bottom for easy mobility to any place.



Model No.	IIC-174	IIC-174 A	IIC-174 B
No. of Shelves	4	5	6
No. of Flourescent Tubes	12	16	20
No. Of Incandescent Lamps	12	16	20

Options:-  
Photoperiodic Timer/Controller to control Cyclic illumination.

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in