




# ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

**IIC-130**



## ROTARY HIGH VACUUM PUMP

- Rotary Oil sealed Pumps are quiet operating & Vibration less.
- MS Angle & sheet housing duly painted.
- Designed for continuous duty.
- Pump mechanism contained in metal case with oil sight gauge drain plug at base.
- Motor Driven Pumps are completely mounted on painted rectangular base plate & furnished with pulley & V belt.

### Applications :

Bio Tech Lab. Pharmaceutical lab. Cosmetic products research lab. Institute.

Model	IIC 130	IIC 130 A	IIC 130 B
Capacity Per/Min	25	50	75
Motor (HP)	1/4 HP	1/4 HP	1/4 HP
Pump Speed	450	450	450

### Technical Parameters :

- Material : MS Angle & sheet plate housing duly powder coated
- Power Consumption : 220 / 230V AC Supply Single Phase 50 HZ
- Duty Cycle : Continuous

### Optional

- Single & Double stage V. Pump (As per Customer requirement capacity of motor & Pump same as above)
- Non Return Valve
- Moisture Trap Vacuum Gauge
- Available in higher capacity also

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in