



ICON INSTRUMENTS COMPANY

IIC-141



VERTEX SHAKER (TEST TUBE SHAKER)

- Vibration type variable speed genie short rapid electronic stroke motor drives.
- Neoprene cup in parts vigorous agitation to tube or small flask.
- Shaking starts when vessel is pressed against cup.
- Fitted with Vacuum Rubber at the bottom to prevent walking.
- Outer casing made of MS /CRC sheet duly to painted / Powder coated

Applications :

For Industrial Research, Bio Tech Medical

Model	IIC 141	IIC 141 A
	Touch On/ Off	Without Touch

Technical Parameters :

- Control Panel : Mains, Mains Indicator, Speed Regulator
 Material : CRC/MS sheet duly Painted
 Power Consumption : 220/230V AC supply Single Phase 50 HZ
 Duty Cycle : Continuous

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in