



ICON INSTRUMENTS COMPANY

(Manufacturer of Laboratory Equipment)

WATER STILL

- Precise pure Pyrogen free distilled water Design feature ensure economy in operation.
- All S. Steel 304 grade body.
- Unique design of water over flow & water level & outlet.
- Kanthal A-1 super quality kettle element for better result / ejection type.
- Condenser tube & Jacket all made of S. Steel 304 grade sheet / pipe.
- Condenser tube can be dismantled by removing the wing nut which holds condenser & Jacket.
- With wall bracket to be fixed on wall on which the still can rest



Applications :

Useful for Hospitals, Medical Research lab. Biotech lab. & Industrial application etc.

Model	IIC 111	IIC 111 A	IIC 111 B
Capacity (L)	2 Liter	4 Liter	6 Liter
Heater Wattage "W"	1.5 K. W.	3 KW. (2 x 1.5Kw.)	4.5 KW, 1.5 KW x 3

Table Model :- All S. Steel Construction as mentioned above Technical parameters remains same

Model II	IIC 111 C	IIC 111 D
Capacity	4 Liter Table Model	6 Liter Table Model

Technical Parameters :

- Duty Cycle : Continuous Check water Level at regular interval
- Power Consumption : 220/230V AC \pm 10%, 50HZ single Phase
- Material : All Stainless Steel 304 grade

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in