



ICON INSTRUMENTS COMPANY

DIGITAL pH METER

- Our Digital PH Meter which is developed using optimum quality raw material like plastic products, PVC, UHM, Brass stainless steel, hastelloy, mild steel and carbon steel. Our Digital pH meters are available in two models and are used to measure acidity or alkalinity of solution. these are widely used in chemical, pharmaceutical, sugar, food cosmetics, plating, water conditioning and waste water treatment industries where pH monitoring is required.

IIC-401



INSTRUMENT TYPE		DIGITAL PH METER
Measurement	pH & ORP	pH & ORP
Display	3 1/2 Digit LED	3 1/2 Digit LED
Range Resolution	0-14.00pH 0 ± 1999mv 0.01pH, 1mV	0-14.00pH 0 ± 1999mv 0.01pH, 1mV
Accuracy	±0.01pH, ±1mV	±0.01pH, ±1mV
Temperature Compensation	Auto + Manual 0 to 100 C	Manual 0 to 100 C
Slope Control	80 to 120 %	80 to 120 %
Recorder Output	0 to 10m V / pH Adjustable	0 to 10m V / pH Adjustable
Power Requirement	220V AC ± 10% 50Hz	220V AC ± 10% 50Hz

DIGITAL CONDUCTIVITY METER

- Leveraging on our high tech manufacturing facilities, we manufacture quality Digital Conductivity Meter. Available at reasonable price and in various specifications, these instruments are widely used by construction industry, steel plants, ready mix concrete plants and cooling towers-chilling, conductivity meters are offered to our clients in standard and customized specifications.
- Conductivity Meter
- Digital Conductivity Mete - Digital LED display with facility for digital Cell, constant adjustment, alongwith Conductivity cell.
- Deluxe TDS Mete - Digital LED display with facility for automatic temperature compensation and Digital cell constant adjustment, along with Conductivity cell and temperature probe.
- Power : 220V +10% 50Hz

IIC-402



Technical Parameters :

INSTRUMENT TYPE	DIGITAL CONDUCTIVITY METER
Display	3 1/2 Digit LED
Range	5 Range 0 to 1000m Mhos/cm
Accuracy	±1% FS± 1 Digit
Temperature Compensation	Manual 0 to 100 C
Cell Constant Adjustment	Adjustable Digital Display
Conductivity / TDS Cell	DIP Type
Power	220V ± 1

Works : Z-67, Okhla Industrial Area, Phase-II, New Delhi - 110020

Phone : +91-11-26383676 • E-mail : iconinstruments.india@gmail.com, sales@iconinstruments.co.in

Website: www.iconinstruments.co.in